



Factors for successful clarification and cell harvesting

Intellectual Property Notice: Part of GE Healthcare's Life Sciences business was acquired by Danaher on 31 March 2020 and now operates under the Cytiva™ brand. Certain collateral materials (such as application notes, scientific posters, and white papers) were created prior to the Danaher acquisition and contain various GE owned trademarks and font designs. In order to maintain the familiarity of those materials for long-serving customers and to preserve the integrity of those scientific documents, those GE owned trademarks and font designs remain in place, it being specifically acknowledged by Danaher and the Cytiva business that GE owns such GE trademarks and font designs.

cytiva.com

GE and the GE Monogram are trademarks of General Electric Company.
Other trademarks listed as being owned by General Electric Company contained in materials that pre-date the Danaher acquisition and relate to products within Cytiva's portfolio are now trademarks of Global Life Sciences Solutions USA LLC or an affiliate doing business as Cytiva.
Cytiva and the Drop logo are trademarks of Global Life Sciences IP Holdco LLC or an affiliate.
All other third-party trademarks are the property of their respective owners.
© 2020 Cytiva
All goods and services are sold subject to the terms and conditions of sale of the supplying company operating within the Cytiva business. A copy of those terms and conditions is available on request. Contact your local Cytiva representative for the most current information.
For local office contact information, visit [cytiva.com/contact](https://www.cytiva.com/contact)

Factors for successful clarification and cell harvesting

This technical brief describes process variables you must consider when selecting a hollow fiber cartridge for clarification and cell harvesting processes.

Four key variables include:

- Membrane pore size
- Hollow fiber path length
- Fiber lumen diameter
- Process conditions

Membrane pore size

Membrane pore size can influence process efficiency in unexpected ways, and specifying a membrane with an appropriate pore size is a key step in developing an effective and efficient clarification or cell harvesting process. For example, when harvesting cells, smaller pore size membranes often provide the highest permeate flux once the system is in a steady state (Fig 1). Hence, experts often select our large pore size ultrafiltration membrane for harvesting *E. coli*, even though it has a relatively small pore size compared to the size of the cells.

Protein passage is paramount in clarification processes, and experienced users typically select open pore size microfiltration membranes (0.2 to 0.65 μm pore size), especially for larger recombinant proteins and monoclonal antibodies. In general, choose a membrane pore size that is at least ten times larger than the target material that you want to pass through the membrane.




Figure 1. Flux of three membranes with all parameters held constant except pore size. Membrane A has a larger pore size than membrane B, which has a larger pore size than membrane C.

In general, cell harvesting is best performed with tight pore size microfiltration membranes or open pore size ultrafiltration membranes (Table 1). For clarification applications, choose a membrane pore size that is at least ten times larger than the target material that you want to pass through the membrane (Table 2).

Hollow fiber path length

The concentration of particulates (cells and cell debris) in the feed material is a primary factor in selecting the length of the cartridge. If the particulate level is high, as is often the case in cell processing applications, you should use short path length cartridges (nominally 30 to 60 cm). Longer path length cartridges (nominally 110 cm) exhibit higher inlet pressure due to increased dynamic friction at a given flow rate.



Fiber lumen diameter

When using hollow fiber cartridges, the inside diameter of the fibers influences other process variables. Feed streams with particulates such as cells and cell debris, flow readily through fibers with larger lumens (inside diameters). For example, hollow fibers with inside diameters of 0.75 to 1.0 mm perform well with particulated feed streams. Likewise, large lumen diameters perform well when processing viscous starting broths and the fouling components often found in upstream processes (for example, lipopolysaccharides in the case of *E. coli*).

Process conditions—cell harvesting

Cell harvesting is typically treated as a dewatering process. During dewatering, broth passes through the hollow fiber membrane and out of the permeate port of the filter. Permeate flows can be high (80-120 l/h) when the feed has a low cell density or lower when the feed has higher concentrations of cells, as with *E. coli*, yeast, or lysate concentration.

Process conditions—clarification

Clarification processes are generally influenced by three key process conditions including recirculation flow rate, permeate flow control, and the timing of diafiltration.

Recommended recirculation flow rate

Experience shows that the following recirculation flow rates work well:

2,000 to 4,000 sec⁻¹ for fragile mammalian cells and viruses
6,000 to 8,000 sec⁻¹ for yeast due to high viscosity
8,000 to 16,000 sec⁻¹ for bacterial cells, lysates, and most proteins

Permeate flow control

Operators must exert deliberate control over transmembrane pressure and careful timing of concentration and cell washing to promote protein passage. During clarification, excessive transmembrane pressure can foul the membrane. Even when retentate pressure is zero, you can reduce transmembrane pressure further by controlling (limiting) the permeate flow from the cartridge. Many users thus add a pump to the permeate line to limit the flow of permeate from the cartridge and thereby exerting some backpressure on the membrane (Fig 2).




Figure 2. Typical cross flow filtration process showing a pump on the permeate line for permeate flow control

How much you reduce the permeate flow (compared to unregulated flow) depends on the nature of the starting feed stream. We recommend lower permeate flow rates: (1) for large target molecules, (2) if the particulates have a broad size distribution, or (3) if the particulates are sticky and fouling (Fig 3). A typical monoclonal antibody clarification from hybridoma cell culture with intact cells may be controlled at approximately 30 to 50 l/h, whereas enzyme clarification from bacterial lysate is almost always at approximately 10 l/h.




Figure 3. Typical permeate flow rates for controlled and uncontrolled process conditions. Permeate flow control results in a more stable flow, higher protein yields, and often shorter process times.

Diafiltration

To promote maximum protein passage during clarification, sound process design is required. For example, particulates can interfere with the passage of protein as the particulates become more concentrated in the feed stream. In addition, fine particles can coat the membrane surface, creating a secondary rejection layer that will retain the target protein. To overcome these problems, a partial concentration should be followed by diafiltration. Proper timing of the diafiltration step is essential. We recommend performing a brief diafiltration at a point where protein is still passing freely—that is, not being retained by the gel layer of concentrated particles that forms on the membrane surface.

Sizing up the process

During the early stages of development, it is important not to “overwhelm” a membrane with a large volume of starting material. You should try to anticipate the relative rate of productivity and introduce a proportionate amount of feed stock to result in a reasonable process time. Here is an example:

When working with a high solids load from yeast fermentation, how much starting material should be used for a one- to two-hour process?

The estimated flux rate is 30 l/mh at room temperature. The membrane surface area is 50 cm² (cartridge model CFP-1-E-2U), and the process design includes a 2× concentration followed by a 3× wash. The calculated permeate flow rate is 150 ml/hr.

If you start with 150 ml and follow the process design, there would be 225 ml of clarified filtrate at the end. Using these estimates, the 150-ml starting volume would be appropriate for a two-hour run.

Summary

There are many factors influencing cross flow filtration clarification. Four key factors include membrane pore size, cartridge length, fiber lumen size, and process conditions.

You must understand these factors to develop an effective and efficient cross flow filtration clarification process.

For details, please contact your GE representative today.

Table 1. Hollow fiber selection guidelines for cell harvesting applications

Factor	Recommendation for cell harvesting
Membrane pore size	500,000 or 750,000 NMWC for ultrafiltration 0.1 or 0.2 µm for microfiltration
Hollow fiber path length	Short path length cartridges (nominally 30 to 60 cm)
Fiber lumen diameter	Large lumen diameter fibers (0.75 to 1.0 mm)
Process conditions	Recirculation flow rate—8,000 to 16,000 sec ⁻¹ shear rate for fouling feed streams, 2,000 to 4,000 sec ⁻¹ for shear-sensitive feed streams. Process sequence—Flux (and protein passage) is dependent on the concentration of particulates. To get the best protein passage, position a cell-washing step before flux declines significantly. Use low to moderate transmembrane pressure (<1 barg [15 psig]). Process temperature—Room temperature works best, but only if process components are stable at this temperature; otherwise 4°C-12°C (39°F-53.6°F) works well, but with lower flux.

Table 2. Hollow fiber selection guidelines for clarification applications

Factor	Recommendation for cell clarification
Membrane pore size	0.2 to 0.65 µm pore sizes for most mammalian cell culture clarifications
Hollow fiber path length	Short path length cartridges (nominally 30 to 60 cm)
Fiber lumen diameter	Large lumen diameter fibers (0.75 to 1.0 mm)
Process conditions	Recirculation flow rate—8,000 to 16,000 sec ⁻¹ shear rate for fouling feed streams, 2,000 to 4,000 sec ⁻¹ for shear-sensitive feed streams. Process sequence—Flux (and protein passage) is dependent on concentration of particulates. To get the best protein passage, position a cell-washing step after a partial concentration, but before flux declines. Use of permeate flow control is strongly recommended. Low transmembrane pressure (<0.7 barg [10 psig]) is usually very important for best protein passage to prevent blinding of membrane by particulates. Process temperature—Room temperature works best, but only if process components are stable at this temperature; otherwise 4°C-12°C (39°F-53.6°F) works well, but with lower flux.

www.gehealthcare.com

Global Headquarters GE Healthcare
Little Chalfont
Buckinghamshire, U.K. HP7 9NA

GE, imagination at work and GE monogram are trademarks of General Electric Company.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. GE Healthcare reserves the right, subject to any regulatory and contractual approval, if required, to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation.

© 2006 General Electric Company - All rights reserved.

GE Healthcare Bio-Sciences AB, a General Electric Company.

GE Healthcare Bio-Sciences AB
Björkgatan 30, 751 84 Uppsala, Sweden

GE Healthcare Europe GmbH
Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare UK Ltd
Amersham Place, Little Chalfont, Buckinghamshire, HP7 9NA, UK

GE Healthcare Bio-Sciences Corp
800 Centennial Avenue, P.O. Box 1327
Piscataway, NJ 08855-1327, USA

GE Healthcare Bio-Sciences KK
Sanken Bldg. 3-25-1, Hyakunincho, Shinjuku-ku,
Tokyo 169-0073, Japan

Asia Pacific Tel +65 6275 1830 Fax +65 6275 1829 Australasia Tel +61 2 9899 0999 Fax +61 2 9899 7511 Austria Tel 01/57606-1619 Fax 01/57606-1627 Belgium Tel 0800 73 888 Fax 02 416 82 06 Canada Tel 800 463 5800 Fax 800 567 1008
Central, East, & South East Europe Tel +43 1 972720 Fax +43 1 97272 2750 Denmark Tel 45 16 2400 Fax 45 16 2424 Finland & Baltics Tel +358 (0)9 512 39 40 Fax +358 (0)9 512 39 439 France Tel 01 69 35 67 00 Fax 01 69 41 96 77
Germany Tel 089 96281 660 Fax 089 96281 620 Greater China Tel +852 2100 6300 Fax +852 2100 6338 Italy Tel 02 27322 1 Fax 02 27302 212 Japan Tel +81 3 5331 9336 Fax +81 3 5331 9370 Latin America Tel +55 11 3933 7300 Fax +55 11 3933 7304
Middle East & Africa Tel +30 210 9600 687 Fax +30 210 9600 693 Netherlands Tel 0800 82 82 82 1 Fax 0800 82 82 82 4 Norway Tel 815 65 555 Fax 815 65 666 Portugal Tel 21 417 7035 Fax 21 417 3184 Russia & other C.I.S. & N.I.S Tel +7 (495) 956 5177
Fax +7 (495) 956 5176 Spain Tel 93 594 49 50 Fax 93 594 49 55 Sweden Tel 018 612 1900 Fax 018 612 1910 Switzerland Tel 0848 8028 12 Fax 0848 8028 13 UK Tel 0800 616928 Fax 0800 616927 USA Tel 800 526 3593 Fax 877 295 8102



imagination at work